

EOGas Cartridge

100% Ethylene Oxide Sterilization

80% less gas than the competition.

With Andersen's match load to cartridge method you never have to worry about tanks of EtO. Andersen Products offers cutting-edge technologies and over 30 years of expertise in sterilization with the only system that offers multiload capability.

► **Match Load Volume:**

Two bag sizes available.

The AN1005 is 18" x 24" bag size with only 5 grams of EtO used.

The AN1006 is 22" x 36" bag size with only 11 grams of EtO used.

► **Cost Effective Kits include:**

25 EOGas 100% EtO cartridges

25 Humidichip® relative humidity stabilizing chips



25 EOGas Dosimeters®

25 corresponding size Sterilization Bags

► **Protection:**

The new diffusion chamber provides extra protection to the operator in case of accidental activation. Its robust design and polyurethane body make the EOGas cartridge one of the most reliable vehicles for releasing the ethylene oxide gas.

Operator safety is insured by the double protection of the trigger via a plastic lock and a safety seal.



The gas goes through 3 barriers before it reaches the venting mechanism.

Double protection trigger.



Two sizes of sterilant cartridges and corresponding bags simplify sterilization. Process loads of various sizes simultaneously. Avoid restrictions of single chamber batch processing.



*leader in Ethylene Oxide
sterilization for over 30 years*



Health Science Park / 3202 Caroline Drive
Haw River, NC 27258 USA
336-376 3000 / Fax 336- 376 8153
www.anpro.com / mailbox@anpro.com